

ECG ELECTRODE



TELECTRODE

Connect with Confidence

Produces ECG electrodes called 'Telectrodes'
which are used for various ECG tests

These are excellent components for sensitive monitoring

Bio Protech Inc.

CONTENT

01

Information

02

02

TELECTRODE

Monitoring ECG Electrodes

03

Diaphoretic Solutions

04

Real-time monitoring during imaging processes

05

Enhanced Electrode Solution for Stress Test

06

03

PROTAB

Diagnostics Resting ECG Electrodes

08

04

PRONEO

Pre-wired Infant Neonatal Electrodes


09

05

Technical Specification

10

The perfect companion for heart health.



With our innovative solutions and products, we can assist healthcare facilities in optimizing their operational costs without compromising on the quality of patient care diagnostics.

Best price, Best performance

Top-notch performance at a low price. An economical choice, but no compromise on performance. With help from Bioprotech, you can meet the challenge by standardizing electrodes without compromising on product performance

Comfortable fit


Emphasizing the design and material selection of the product, the ECG electrodes can touch the skin smoothly and be worn comfortably.

Accurate and reliable measurement

High-performance sensor and ECG electrodes of proven quality provide accurate ECG measurement. Quick and consistently excellent signal quality is obtained.

Backing choice

You can choose from 3 types(Foam, Cloth, Clear tape) of backing for conformability, patient comfort and low profile.

	PE Foam The foam material provides a soft and cushioning effect that gently touches the skin. It prevents unnecessary moisture build-up on the skin and supports skin health.	T715	T716	T717	T718	T712T	T711T	T711CT	T724T	T724CT
	Cloth The cloth material adheres reliably to the skin and remains in place during movement. The electrodes provide the long term comfort and durability required by active patients.	T815	T815C	T816	T832	T832C				
	Clear Tape The clear tape material is transparent, and when the ECG electrodes are attached to the skin, the condition under the skin can be visually checked. This can be helpful in monitoring skin conditions or monitoring skin changes.	T916	T916C							

By choosing Bio Protech as your partner, you can benefit from our expertise, cost-effective solutions, and commitment to high-quality patient care diagnostics. Let us help you navigate the challenges of reducing operating costs while maintaining excellence in healthcare.

TELECTRODE

Monitoring ECG Electrodes

Multiple sizes and shapes

We offer ECG electrodes for adult patients in a variety of sizes and shapes to suit different body shapes and sizes.

ECG Electrodes for Universal

Model	Backing	Size(mm)	Shape
T715	PE Foam	Dia. 55	Round
T716	PE Foam	43 x 45	Teardrop
T717	PE Foam	Dia. 43	Round
T711T	PE Foam	37 x 42.5	Teardrop
T712T	PE Foam	45 x 49	Teardrop
T724T	PE Foam	36 x 42	Square
T816	Cloth	43 x 45	Teardrop
T916	Clear Tape	43 x 45	Teardrop



T716



T724T



T916

Small size option

The T718 electrode is compact size and design to suit the pediatric's anatomy.

ECG Electrodes for Pediatric

Model	Backing	Size(mm)	Shape
T718	PE Foam	32 x 38	Teardrop



T718

Suggested Applications	T715	T716	T717	T718	T711T	T712T	T724T	T816	T916
EMS/ED	○	○	○	○				○	○
ICU/CCU	○	○	○	○				○	○
Telemetry/Holter/Event					○	○	○	○	
OR	○	○	○	○	○	○	○		○
Cath Lab									
Stress Test							○		
Cardiac Rehab	○	○	○	○	○	○	○	○	
Echo-Cardio									
Universal Foam Monitoring	○	○	○	○			○	○	
Pediatric				○					
Solutions of Fragile skin									○
Diaphoretic solutions									
Fluid resistant solutions	○	○	○	○		○	○		○

TELECTRODE

Diaphoretic Solutions

Expansive adhesive surface for enhanced contact and adherence.

TELECTRODE T815 electrodes integrate an expansive surface area and a robust adhesive to enhance the quality of trace measurements.

Model	Backing	Size(mm)	Shape
T815	Cloth	Dia. 55	Round



T815

Solutions for Fragile Skin

Soft adhesive

Diaphoretic Monitoring Electrodes enhance the quality of the trace by offering a larger surface area for better contact and a powerful adhesive to ensure secure attachment.

Skin-friendly preparatory abrader for a gentle skin surface.

The use of a gentle abrader in TELECTRODE ensures a soft and smooth skin surface for optimal EKG trace quality.



T832

Model	Backing	Size(mm)	Shape
T832	Cloth	Dia. 32	Round
T832C	Cloth	Dia. 32	Round

Suggested Applications	T815	T832	T832C
EMS/ED	○	○	○
ICU/CCU	○	○	○
Telemetry/Holter/Event			
OR			
Cath Lab			○
Stress Test			
Cardiac Rehab	○	○	○
Echo-Cardio			○
Universal Foam Monitoring			
Pediatric		○	
Solutions of Fragile skin		○	○
Diaphoretic solutions	○		
Repositionable		○	○

TELECTRODE

Real-time monitoring during imaging processes

Saves Time and Cost.

With the electrodes already applied in other departments, there is no need to remove them during X-ray or MRI procedures. This eliminates the hassle of changing electrodes and streamlines the process for both healthcare professionals and patients.

Model	Backing	Size(mm)	Shape
T715C	PE Foam	Dia. 55	Round
T716C	PE Foam	43 x 45	Teardrop
T718C	PE Foam	32 X 38	Teardrop
T815C	Cloth	Dia. 55	Round
T832C	Cloth	Dia. 32	Round
T711CT	PE Foam	37 X 42.5	Teardrop
T724CT	PE Foam	36 X 42	Square
T916C	Clear Tape	43 x 45	Teardrop



T718C



T724CT

Suggested Applications	T715C	T716C	T718C	T815C	T832C	T711CT	T724CT	T916C
EMS/ED	○	○	○	○	○			○
ICU/CCU	○	○	○	○	○			○
Telemetry/Holter/Event						○	○	
OR	○	○	○			○	○	○
Cath Lab	○	○	○	○	○	○	○	○
Stress Test						○	○	
Cardiac Rehab	○	○	○	○	○	○	○	○
Echo-Cardio	○	○	○	○	○	○	○	○
Universal Foam Monitoring							○	
Pediatric			○		○			
Solutions of Fragile skin							○	○
Diaphoretic solutions								
Fluid resistant solutions	○	○	○			○	○	○

TELECTRODE

Enhanced Electrode Solution for Stress Test and Holter

Strong signals and Stable baselines.

The conductive wet gel rapidly and effectively reduces skin impedance.


The wet gel enhances the contact between the skin and electrode.

Utilizing both quick and durable adhesives, the electrode ensures reliable fixation and signal stability.


The Ag/AgCl sensor guarantees consistent signal quality.

ECG Electrodes for Stress Test and Holter

Model	Backing	Size(mm)	Shape
T715W	PE Foam	Dia. 55	Round
T716W	PE Foam	43 x 45	Teardrop
T815W	Cloth	Dia. 55	Round
T816W	Cloth	43 x 45	Teardrop



T716W



T815W

Suggested Applications	T715W	T716W	T815W	T816W
EMS/ED				
ICU/CCU				
Telemetry/Holter/Event	○	○	○	○
OR				○
Cath Lab				
Stress Test	○		○	
Cardiac Rehab				○
Echo-Cardio				
Universal Foam Monitoring				
Pediatric				
Solutions of Fragile skin				
Diaphoretic solutions				
Fluid resistant solutions				

TELECTRODE

Enhanced Electrode Solution for Stress Test and Holter

Strong signals and Stable baselines & Radiolucent

The conductive wet gel rapidly and effectively reduces skin impedance.

The wet gel enhances the contact between the skin and electrode.

Utilizing both quick and durable adhesives, the electrode ensures reliable fixation and signal stability.

The Ag/AgCl sensor guarantees consistent signal quality.



T715CW



T815CW

ECG Electrodes for Stress Test and Holter & Carbon

Model	Backing	Size(mm)	Shape
T715CW	PE Foam	Dia. 55	Round
T716CW	PE Foam	43 x 45	Teardrop
T815CW	Cloth	Dia. 55	Round
T816CW	Cloth	43 x 45	Teardrop

Suggested Applications	T715CW	T716CW	T815CW	T816CW
EMS/ED				
ICU/CCU				
Telemetry/Holter/Event	○	○	○	○
OR				
Cath Lab	○	○	○	○
Stress Test	○		○	
Cardiac Rehab				
Echo-Cardio				
Universal Foam Monitoring				
Pediatric				
Solutions of Fragile skin				
Diaphoretic solutions				
Fluid resistant solutions				

PROTAB

Diagnostics Resting ECG Electrodes

The unique needs of resting ECG monitoring are met with the Bioprotech family of electrodes designed for this critical diagnostic application. Whether your preference is for a tab or snap connection, Bioprotech can fulfill the need.

Easy to apply and remove.

With this design, accessing the electrode becomes quick and effortless. The tab itself is strategically positioned for convenient leadwire placement and ensures easy removal when needed.

Intimate skin contact.


The conductive adhesive used in PROTAB Electrodes establishes a close and precise connection with the skin, guaranteeing accurate EKG traces, even in difficult-to-reach areas.

Superior adhesion.

The electrode's excellent adhesion allows for secure placement, while still providing the flexibility to reposition it as needed. When a higher level of adhesion is required, PROTAB is the preferred choice.

Defibrillation recoverable.

The PROTAB electrodes can withstand the effects of defibrillation and remain intact, allowing for continued use without the need for replacement.



PT1334 |
13x34(mm)
10pcs/1card
Radiolucent



PT2034 |
20x34(mm)
10pcs/1card
Radiolucent



PT2334 |
23x34(mm)
10pcs/1card
Radiolucent



PT3034 |
30x34(mm)
10pcs/1card
Radiolucent

Suggested Applications	PT1334	PT2034	PT2334	PT3034
EKG Lab	☐	☐	☐	☐
12 Lead EKG	☐	☐	☐	☐
Short Term Diagnostics	☐	☐	☐	☐
Pediatric	☐	☐		



PRONEO

Pre-wired Infant Neonatal Electrodes

Bioprotech offers a specialized range of pre-wired electrodes tailored specifically for newborns and infants. The PRONEO is meticulously designed to meet the unique needs of the youngest patients, ensuring optimal comfort and accuracy during monitoring.

Pre-wired and color-coded.


For ease of application and trace quality.

Gentle adhesive.

Minimizes the potential for skin injury on developing neonatal skin prior to complete barrier formation.


Smaller footprint for neonates.

Accommodates smaller chest sizes.



SG15 | Socket $\Phi 1.5$
Electrode 22mm
Y.G.R
Lead-wire 60cm

Suggested Applications	SR15	SR20	SR15L	SR20L	SR15S	SR15Limb
Echo-Cardio	○	○	○	○	○	○
Peds ICU			○	○		
NICU	○	○			○	○
Peds OR			○	○		
Pediatric			○	○		
Prewired Infant & Neonatal	○	○			○	○



SR20 | Socket $\Phi 2.0$
Electrode 22mm
W.B.G
Lead-wire 60cm
Radiolucent

Suggested Applications	SG15	SG20	SG15L	SG20L	SG15S	SG15Limb
Echo-Cardio	○	○	○	○	○	○
Peds ICU			○	○		
NICU	○	○			○	○
Peds OR			○	○		
Pediatric			○	○		
Prewired Infant & Neonatal	○	○			○	○



SR15L | Socket $\Phi 2.0$
Electrode 22x36mm
White
Lead-wire 60cm
Radiolucent



Technical Specification

TELECTRODE

If you need a strip type, consult with our sales team.

Model	Backing	Size(mm)	Shape	Description
T715	PE Foam	Ø 55	Round	Adult
T715C	PE Foam	Ø 55	Round	Adult
T716	PE Foam	43 x 45	Teardrop	Universal
T716C	PE Foam	43 x 45	Teardrop	Universal/ Radiolucent
T717	PE Foam	Ø 43	Round	Adult
T718	PE Foam	32 x 38	Teardrop	Pediatric
T718C	PE Foam	32 x 38	Teardrop	Pediatric / Radiolucent
T815	Cloth	Ø 55	Round	Adult
T815C	Cloth	Ø 55	Round	Adult / Radiolucent
T816	Cloth	43 x 45	Teardrop	Universal
T916	Claer tape	43 x 45	Teardrop	Universal
T916C	Claer tape	43 x 45	Teardrop	Universal/ Radiolucent
T832	Cloth	Ø 32	Round	Pediatric
T832C	Cloth	Ø 32	Round	Pediatric / Radiolucent
T711T	PE Foam	37 x 42.5	Teardrop	Universal
T711CT	PE Foam	37 x 42.5	Teardrop	Universal / Radiolucent
T712T	PE Foam	45 x 49	Teardrop	Universal
T724T	PE Foam	36 x 42	Square	Universal
T724CT	PE Foam	36 x 42	Square	Universal / Radiolucent
T715W	PE Foam	Ø 55	Round	Adult / Wet gel
T715CW	PE Foam	Ø 55	Round	Adult / Wet gel / Radiolucent
T716W	PE Foam	43 x 45	Teardrop	Universal / Wet gel
T716CW	PE Foam	43 x 45	Teardrop	Universal / Wet gel / Radiolucent
T815W	Cloth	Ø 55	Round	Adult / Wet gel
T815CW	Cloth	Ø 55	Round	Adult / Wet gel / Radiolucent
T816W	Cloth	43 x 45	Teardrop	Universal / Wet gel
T816CW	Cloth	43 x 45	Teardrop	Universal / Wet gel / Radiolucent

PACKING | 50pcs/pack, 20packs/box, 10boxes/case

WETGEL PACKING | 30pcs/pack, 20packs/box, 10boxes/case

PROTAB


Model	Shape	Size(mm)	Material	Patient	Use
PT1334	Poise	13 x 34 (mm)	PET, Ag/AgCl	General	Diagnostic
PT2034	Poise	20 x 34 (mm)	PET, Ag/AgCl	General	Diagnostic
PT2334	Poise	23 x 34 (mm)	PET, Ag/AgCl	General	Diagnostic
PT3034	Poise	30 x 34 (mm)	PET, Ag/AgCl	General	Diagnostic

PACKING | 10pcs/card, 10cards/pack, 5packs/box, 5boxes/case

PRONEO

Model	Socket size	Electrode size(mm)	Radiolucent	Wire length	Backing
SG15	Φ1.5	22 x 22	x	60 cm leadwire	Cloth
SG20	Φ2.0	22 x 22	x	60 cm leadwire	Cloth
SG15L	Φ1.5	22 x 36	x	60 cm leadwire	Cloth
SG20L	Φ2.0	22 x 36	x	60 cm leadwire	Cloth
SG15S	Φ1.5	15 x 20	x	60 cm leadwire	Cloth
SG15Limb	Φ1.5	17 x 22	x	60 cm leadwire	Cloth
SR15	Φ1.5	22 x 22	O	60 cm leadwire	Cloth
SR20	Φ2.0	22 x 22	O	60 cm leadwire	Cloth
SR15L	Φ1.5	22 x 36	O	60 cm leadwire	Cloth
SR20L	Φ2.0	22 x 36	O	60 cm leadwire	Cloth
SR15S	Φ1.5	15 x 20	O	60 cm leadwire	Cloth
SR15Limb	Φ1.5	17 x 22	O	60 cm leadwire	Cloth

PACKING | 3pcs(set)/pack, 100packs/box, 12boxes/case



Sales office

Address | Orange County Room 1010~1011, 993, Seongnam-daero, Jungwon-gu, Seongnam-si, Gyeonggi-do, 13438 Republic of Korea

Mail | info@protechsite.com **Phone** | +82-33-735-7720 **Fax** | +82-70-4896-5483

About Us

Head Office | 151-3, Donghwagongdan-ro, Munmak-eup, Wonju-si, Gangwon-do, 26365 Korea

E-mail | info@protechsite.com **Phone** | +82-33-735-7720 **Fax** | +82-33-735-7736

North American Headquarter | 13311-a 166th St. Cerritos, CA 90703, USA

Phone | 909-464-1464 **Fax** | 909-464-2664)



BIO PROTECH